

Welcome to SeqBox

Mini-tutorial to get started
full tutorial at:

<https://kendomaniac.github.io/docker4seq/index.html>


For support please contact: raffaele.calogero@unito.it


SeqBox login


- Local ethernet IP: 192.169.1.4
 - User: user
 - Password: 4SeqGUIv2
-
- Scratch folder: /data/scratch
 - Genomes folder: /data/genomes
 - In genomes are already available: hg38star, mm10star and hg38salmon

Login via VNC


- Install VNC client
- Connect with an ethernet cable your PC with SeqBox:
- Open VNC
- See next slide









Step3: press return in case you get this screen saver



Step4: type password



Step4: start 4SeqGUI



For help see: <https://kendomaniac.github.io/docker4seq/articles/docker4seq.html>